
**Operating Instructions | Betriebsanleitung |
Installation et mode d'emploi**

Sartorius YDP10-OCE

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English – page 3

In cases involving questions of interpretation,
the German-language version shall prevail.

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Contents

Intended Use	4	Troubleshooting Guide	9
Safety Instructions	4	Maintenance	10
Installation	5	Accessories	14
Operation	8	Specifications	14
Printing a Measured Value from the Connected Weighing Instrument	8	Electromagnetic Compatibility	15
Configuration	8	Using the Printer in Legal Metrology	15
		Test Certificate	16

Intended Use

Model YDP 10-OCE is a data printer designed to receive weighing data from a connected weighing instrument and print it on plain paper (roll width: 57.5 mm). Text and data are printed with up to 24 characters per line.

The data printer meets the quality standards for research in the chemical and pharmaceutical industries, in compliance with the requirements of GLP (Good Laboratory Practice) and GMP (Good Manufacturing Practice).

The printer can be used for output of weighing data. Data is sent to the output port either automatically or at the press of a key. The printer can be used as auxiliary equipment with a weighing instrument that has been verified for use in legal metrology; see “Using the Printer in Legal Metrology” for more information.

Safety Instructions

The data printer complies with the EC Directives and Standards for electrical equipment and electromagnetic compatibility.

Use only rechargeable nickel-metal hydride, size AA (order no. 69Y03136). Avoid draining the batteries completely, as this will reduce their recharging capacity. Used batteries are hazardous waste and must be properly disposed of in accordance with your local hazardous waste disposal regulations.

As a rule, the printer is powered over the data cable connected to the weighing instrument.

Depending on your weighing instrument model, it may be necessary to use a different external power source. To do this, use only an original Sartorius AC adapter (see “Accessories”). If you have any difficulties, contact our Service Center or your Sartorius dealer.

Use only a moist cloth to wipe down the printer to ensure that no liquid enters the printer housing during cleaning.

Caution: The serrated edge for tearing off printouts is sharp; direct contact can result in injury.

Installation



Equipment Supplied

The following components are supplied:


- Interface cable for connecting the printer to a weighing instrument
- Ink ribbon
- 4 rechargeable batteries
- Paper roll

Installing the Rechargeable Batteries

- Remove the 4 batteries from the packaging (they are discharged)
- Lift the cover of the paper roll compartment
- To open the battery compartment, tilt back the cover.

- Install the 4 rechargeable batteries in the battery compartment. Make sure that the poles match the + and - signs.
- Close the battery compartment
- Close the paper roll compartment
- You cannot thread the paper roll until the printer has been turned on.


Connecting the Data Printer to a Weighing Instrument

- Connect the interface cable to the printer and to the weighing instrument. Secure the connection.

Turning On the Weighing Instrument and the Printer

- Plug in the AC adapter to power the weighing instrument
- Turn on the weighing instrument
 - > The printer switches on automatically

- > The green LED on the left-hand side of the front panel lights up


Charging the Batteries during Initial Operation

- Charge the batteries by leaving the weighing instrument turned on
 - > The batteries are charged while the weighing instrument is in standby mode
 - > It takes 28 hours to charge the batteries to full capacity
- The LED flashes red to indicate that there is not enough charge in the batteries to power the printer. If the LED does not stop flashing after a maximum of 3 hours charging, the batteries must be replaced (see page 12).

Installing the Ink Ribbon Cassette


- For details on installing the ink ribbon cassette, see the section entitled "Maintenance" (page 11).
- Inserting the Paper Roll

- For details on inserting the paper roll, see the section entitled “Maintenance” (page 10).

- Plug the original Sartorius AC adapter into an electrical outlet
- Insert the right-angle plug into the jack on the right-hand side of the printer's rear panel. The green LED lights up to indicate the presence of voltage.
- For details on installing the batteries, see page 6
- Follow the steps described in the previous section to connect the printer and the weighing instrument

Operation

The data printer has 2 function keys for operation:

- (): Line feed; advances the paper by one blank line.
- (): Print measured value or begin automatic printing

Printing a Measured Value from the Connected Weighing Instrument

- Advance paper to desired position, if necessary:
Press the () key (repeatedly if desired), or press and hold for continuous line feed
 - To print data, press the () key
- > A sample printout is shown here on the left

| 001: + 123.4 g |

Configuration

Configuring the Operating Menu

Note:

The connected weighing instrument must be configured for data communication with the printer.



Configuring the Weighing Instrument's Printer Port

The following settings are configured in the printer at the factory:

9600 baud, 8 bits, odd parity, one stop bit, hardware handshake (CTS/DTR)

When the printer is turned on, these settings are loaded in the printer's memory (EEPROM) for use during operation.

Troubleshooting Guide


LED	Cause	Solution
Red	Data in printer memory (EEPROM) defective (checksum invalid)	<ul style="list-style-type: none">● Press the () key> Factory settings are loaded○ It may be necessary to configure menu settings for the printer port and in the statistics application
Red	The capacity of the input buffer has been exceeded	<ul style="list-style-type: none">● Switch the weighing instrument off and then on again● Press the () key
Red	No control signal from printer	<ul style="list-style-type: none">● Press the () key● Check the batteries○ Recharge batteries if necessary○ Replace the batteries if necessary
Flashing red	Low battery power Batteries dead	<ul style="list-style-type: none">● Recharge batteries (see “Maintenance”)
Red	Transmission error in a data string	<ul style="list-style-type: none">● Make sure the interface parameter settings in the printer match those in the weighing instrument● Press the () key

Maintenance

Replacing the Paper Roll

The paper rolls used with your data printer must have the following specifications (see “Accessories”):


- Standard paper roll
- Wood-free paper, 60 g
- Width: 57.5 ± 0.5 mm
- Diameter: Max. 65 mm
- Core diameter: Min. 8 mm


- Lift the cover of the paper roll compartment
- Pull out the roll holders from both sides of the roll core as shown on the left

- Remove the old roll core

- Place a new paper roll on a flat surface
- The open end must be on the bottom of the roll and facing the front of the printer


- Hold the paper roll above the paper roll compartment between thumb and index finger of one hand
- Thread the paper through the slot underneath the ink ribbon with your other hand
- To transport the paper through the ink ribbon cassette, press the () key repeatedly, or hold it down
- > The paper will be threaded into, then out of the ink ribbon cassette

- To align the paper roll, hold it so that you can push both roll holders into the core

- Release the roll holders on both sides
- Close the paper roll compartment


Replacing the Ink Ribbon Cassette
Your data printer requires a black ink ribbon cassette, type ERC-22 (see “Accessories”).

- Lift the cover of the paper roll compartment
- To remove the paper from the ink ribbon cassette, turn the paper roll opposite to the paper feed direction. Press down on the right-hand side of the cassette.

> The left-hand side of the cassette rises

- Remove the cassette


- Insert a new ink ribbon cassette from the top

> The cassette locks into place

- Turn the small tension wheel located on the left in the direction of the arrow until the ribbon is stretched taut in front of the paper


- To thread the paper, see “Replacing the Paper Roll”
- Close the paper roll compartment

Charging the Batteries

If the LED is red (not flashing), the batteries are completely drained.

- The LED flashes red to indicate that the batteries are charging; do not switch off the weighing instrument.
- > The batteries are charged while the weighing instrument is in standby mode
- > It takes 28 hours to charge the batteries to full capacity
- If the LED does not stop flashing red within a maximum of 3 hours, replace the batteries

Replacing the Batteries

⚠ Avoid draining the batteries completely, as this will reduce their recharging capacity.


⚠ Use only size AA rechargeable nickel-metal hydride batteries (see "Accessories"; order no. 69Y03136).

- Lift the cover of the paper roll compartment
- To open the battery compartment, tilt back the cover

⚠ Remove the 4 batteries

Do not dispose of batteries in normal household waste. Rechargeable batteries contain toxic materials and must be disposed of in accordance with local waste disposal regulations.

- If you have any questions, please contact your local Sartorius Service Center

- Place the 4 new batteries in the battery compartment. Make sure that the poles match the + and - signs
- Press down on the cover to close the battery compartment
- Close the paper roll compartment

Cleaning the Data Printer

- ⚠ When cleaning the data printer, make sure that no liquid enters the housing
- ⚠ Do not use any aggressive cleaning agents (solvents or similar agents).
- Disconnect the power supply by unplugging the interface cable (connecting the weighing instrument) from the rear panel of the printer.
- If you are using an external power supply, unplug the AC adapter from the electrical outlet (mains supply)
- Use only a cloth moistened with a mild detergent (soap) to wipe down the printer housing
- After cleaning, wipe down the printer with a soft, dry cloth

Accessories

Product	Order No.
Spare paper rolls (box of 5)	6906937
Ink ribbon cassette	6906918
AC adapter with specifications rated for:	
Europe	6971412
UK	6971414
USA	6971413
Australia	6971411
South Africa	6971410
Set of 4 rechargeable batteries (nickel-metal hydride, size AA)	69Y03136
Data cable for connecting the printer to a weighing instrument or other device of an older design series (MP technology), for use in residential areas (EN55022, Class B)	VF1503 (Fast Factory)

Specifications

Dimension	Value
Length	230 mm
Width	158 mm
Height	80 mm
Weight (without batteries or paper roll)	Approximately 800 g
Operating temperature	0°C to +40°C (32°F to 104°F)
Allowable storage temperature	-40°C to +70°C (-40°F to +158°F)

Electromagnetic Compatibility

CE Marking

The printer complies with the following EC Directives and European Standards: 2004/108/EC “Electromagnetic Compatibility (EMC)” Applicable European Standards: Limitation of emissions: In accordance with product standard EN 61326-1 Class B (residential area) Defined immunity to interference: in accordance with product standard EN 61326-1 (industrial environment)

Note: The operator shall be responsible for any modifications to Sartorius equipment and for any connections of cables or equipment not supplied by Sartorius and must check and, if necessary, correct these modifications and connections. On request, Sartorius will provide information on the minimum operating specifications (in accordance with the standards listed above for defined immunity to interference).

Council Directive 2006/95/EC: “Electrical equipment designed for use within certain voltage limits”

Applicable European Standards: EN 61010-1 Safety requirements for electrical equipment for measurement, control and laboratory use Part 1: General requirements If you use electrical equipment in installations and under ambient conditions subject to stricter safety standards than those described in the manual, you must comply with the provisions as specified in the applicable regulations for installation in your country.

Using the Printer in Legal Metrology

The German National Institute of Physics and Metrology (PTB) has issued an independent test certificate for the YDP10-OCE (see certificate for number). With this certificate, the printer is permitted for use in legal metrology in the EU and the EEA. This applies only when the printer is connected to a weighing instrument that has an EC type approval certificate, or a national equivalent, verifying that the weighing instrument is verifiable for legal metrology.

If the weighing instrument has one or more display digits marked ($e \neq d$), the printer must mark these digits as well; for example, with brackets or inverse printing. This does not apply for weighing instruments of accuracy class I with a verification scale interval $e = 1 \text{ mg}$ and a scale interval of $d < 0.1 \text{ mg}$. For weighing instruments with $e = d$, marking digits is not permitted.

Data output from the weighing instrument must be configured so that the marking of digits is passed to the printer and included on the printout.

If the printer is subsequently connected to a verified weighing instrument already in use in legal metrology, notify your local weights and measures inspection office or contact your local Sartorius service center if you wish to have the data printer configured and prepared for verification for use in legal metrology. The above does not concern printers that are immediately used in legal metrology.

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Object:

serieller Schnittstelle

Modul / Printer to be connected to verified non-automatic weighing

instruments with serial interface

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Copyright by Sartorius AG,
Goettingen, Germany.
All rights reserved. No part
of this publication may
be reprinted or translated in
any form or by any means
without the prior written
permission of Sartorius AG.
The status of the information,
specifications and illustrations
in this manual is indicated
by the date given below.
Sartorius AG reserves the
right to make changes to the
technology, features,
specifications and design of the
equipment without notice.

Status:
August 2008, Sartorius AG,
Goettingen, Germany